INVACARE Yes, you can:

Tel: 01656 776222 Fax: 01656 776220 Email: ordersuk@invacare.com Online Spares available at: www.invacare.co.uk	APRIL 2020 Customer Ref:
Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Kuschall K-Series 2.0

Standard Features

- Max User Weight 130kg (20.5 stone) .
- Crash Tested ISO 7176-Part 19
- Active Rigid Wheelchair with Height Adjustable Rigid • Footrest
- Multiple adjustments, Front and Rear Seat Heights and Centre of Gravity
- Full Hydroformed Aluminium Front Frame 75° Front • Angle
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims and Schwalbe Black and Grey Pneumatic Tyres
- Choice of 10 Seat Widths 320mm to 500mm •
- Choice of 8 Seat Depths 350mm to 525mm •
- Range of Front Seat Heights from 450mm to 520mm •
- Range of Rear Seat Heights from 390mm to 490mm •
- Lockdown Folding, Angle and Height Adjustable •
- Backrest
- Black Nylon with Reflective Trim, Tension Adjustable . Upholstery
- Without Push Handles, Composite Clothesguards •
- Black Aluminium Axle Tube with 3° Camber
- Black Comfort Front Castor .
- Standard Brake Lever Bent
- Tool Kit

Standard Base price: £2,095



Please make your selection with a	cross in the box	0	Free of charge	•	standard feature POA special request
Customer weight =	KG				

Code	Description	
BASE PRODUCT*		
DKA0067	KÜSCHALL K-S	
DKA0000	MYKÜSCHALL	

BASE PRODU		
DKA0067	KÜSCHALL K-SERIES	£2,095
DKA0000	MYKÜSCHALL Special Request - Please include quotation reference in special description at the end of the form where applicable	POA

FRAME*			
DKA1040	V Style Standard Frame Angle 75° Aluminium	(NC 1475, 1490, 1495)	•
DKA1050	Tapered Style Standard Frame Angle 75° Aluminium	(NC 1475, 1490, 1495, 1740, 1745, 1750, 3230)	0
DKA1060	V Style +5cm Frame Angle 75° Aluminium. Frame to accommodate a longer seat depth and seat height for taller clients	(NC 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760)	0
DKA1070	V Style Standard Frame Angle 90° Aluminium	(NC 1475, 1490, 1495, 5300, 5350, 5400, 5450, 5500)	0
DKA1080	Tapered Style Standard Frame Angle 90° Aluminium. Taper before the frame bend creating a parallel style front frame with a wider footrest space	(NC 1475, 1490, 1495, 1740, 1745, 1750, 3230, 5300, 5350, 5400, 5450, 5500)	0
DKA1090	V Style +5cm Frame Angle 90° Aluminium. Frame to accommodate a longer seat depth and seat height for taller clients	(NC 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760, 5300, 5350, 5400, 5450, 5500)	O
SUSPENSION			
DKA1200	Seat Suspension, Weight Range 40-65 KG	(NC 1310, 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430)	£359
DKA1210	Seat Suspension, Weight Range 65-85 KG	(NC 1310, 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430)	£359
DKA1220	Seat Suspension, Weight Range 85-100 KG	(NC 1310, 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430)	£359
SEAT WIDTH*			
DKA1310	320 mm	(NC 2050, 2610, 2630, 2720, 2730)	0
DKA1315	340 mm	(NC 2720, 2730)	0
DKA1320	360 mm		0
DKA1325	380 mm		0
DKA1330	400 mm		•
DKA1335	420 mm		0
DKA1340	440 mm	(NC 2630)	0
DKA1345	460 mm	(NC 2630)	0
DKA1350	480 mm	(NC 2630, 3230)	0
DKA1355	500 mm	(NC 2630, 3210, 3230)	0
SEAT DEPTH*			
DKA1410	350 mm	(NC 1500, 1510, 2630, 2720, 2730)	0
DKA1420	375 mm	(NC 2720, 2730)	0
DKA1435	400 mm		•
DKA1450	425 mm		0
DKA1465	450 mm	(NC 2630)	0
DKA1475	475 mm	(NC 2630)	0
DKA1490	500 mm	(NC 2630)	0
DKA1495	525 mm	(NC 2630, 2610)	0
SEAT POSITION	I. Seat module can be moved forward and backward along	g the frame to achieve a more active or stable rear whee	el position*
DKA1500	Position 1 = 140 mm MOST ACTIVE POSITION		0
DKA1510	Position 2 = 125 mm		0
DKA1520	Position 3 = 110 mm		•
DKA1530	Position 4 = 95 mm		0
DKA1540	Position 5 = 80 mm LEAST ACTIVE POSITION		0

BACK HEIGHT.	Transportation kit is only possible with backrest heights	s of 405mm and above*	
DKA1600	Back Height for 3rd party backrest. Final height depends on the height of the positioning backrest, select if Matrx/Fixed backrest is required.	(NC 2100, 2220, 2230, 7600)	0
DKA1605	270 mm	(NC 2150, 2240, 7600)	0
DKA1610	285 mm	(NC 2150, 2240, 7600)	0
DKA1615	300 mm	(NC 2150, 2240, 7600)	0
DKA1620	315 mm	(NC 2150, 2240, 7600)	0
DKA1625	330 mm	(NC 2150, 2240, 7600)	0
DKA1630	345 mm	(NC 2150, 2240, 7600)	0
DKA1635	360 mm	(NC 2150, 2240, 7600)	0
DKA1640	375 mm	(NC 2150, 2240, 7600)	О
DKA1645	390 mm	(NC 2150, 2240, 7600)	0
DKA1650	405 mm	(NC 2150, 2240)	0
DKA1655	420 mm	(NC 2150, 2240)	0
DKA1660	435 mm	(NC 2150, 2240)	0
DKA1665	450 mm	(NC 2150, 2240)	0
DKA1670	465 mm	(NC 2150, 2240)	0
DKA1675	480 mm	(NC 2150, 2240)	0
FRONT SEAT H	EIGHT. For 75° frames, the maximum lower leg must be	front seat height minus 30mm and minus 40mm for 90° frame	S*
DKA1740	450 mm	(NC 5300, 5350, 5400, 5450, 5500)	О
DKA1745	460 mm	(NC 5300, 5350, 5400, 5450, 5500)	О
DKA1750	470 mm	(NC 5300, 5350, 5400, 5450, 5500)	0
DKA1755	480 mm		0
DKA1760	490 mm		0
DKA1765	500 mm		٠
DKA1770	510 mm		0
DKA1775	520 mm		0
REAR SEAT HE	IGHT*		
DKA1815	390 mm	(NC 3490, 3495, 6030, 6040, 6430)	0
DKA1820	400 mm	(NC 3490, 3495, 6030, 6040)	0
DKA1825	410 mm	(NC 3490, 3495)	0
DKA1830	420 mm		0
DKA1835	430 mm		0
DKA1840	440 mm		٠
DKA1845	450 mm		0
DKA1850	460 mm		0
DKA1855	470 mm		0
DKA1860	480 mm		0
DKA1865	490 mm	(NC 6430)	0

LOWER LEG LI	ENGTH. Cushion height must be deducted*		
DKA1918	290 mm		О
DKA1920	300 mm		O
DKA1922	310 mm		O
DKA1924	320 mm		C
DKA1926	330 mm		O
DKA1928	340 mm		O
DKA1930	350 mm		C
DKA1932	360 mm		O
DKA1934	370 mm		0
DKA1936	380 mm		O
DKA1938	390 mm		O
DKA1940	400 mm		O
DKA1942	410 mm		O
DKA1944	420 mm		O
DKA1946	430 mm		•
DKA1948	440 mm		O
DKA1950	450 mm		O
DKA1952	460 mm		O
DKA1954	470 mm		O
DKA1956	480 mm		O
	TERIAL*		
DKA2000	Aluminium		•
DKA2050	Titanium	(NC 2150, 2430, 7600)	£170
BACKREST TY	PE*		
DKA2100	Standard Backrest with Backbrace	(NC 2240, 2440)	•
DKA2150	Backrest With Lower Backbrace. Backrest Tube Height is 285 mm	(NC 2220, 2230, 2430)	0
BACKREST UP			
DKA2220	Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener	(NC 2440)	•
DKA2230	Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener	(NC 2440)	£185
DKA2240	Without Backrest Upholstery Select if fitting a Matrx/Fixed backrest, no upholstery required		0
PUSH HANDLE	S*		
DKA2400	Without Push Handles		•
DKA2410	Push Handles Fixed Short/Mini Hand Grip		£110
DKA2420	Push Handles Fixed Long Hand Grip		£110
DKA2430	Push Handles Height Adjustable, Rear Set		£210
DKA2440	Push Handles Height Adjustable, Integrated in backrest		£181
DKA2450	Push Handles Foldable		£164

BACKREST AN	GLE*		
DKA2500	74° angle from the seat		0
DKA2520	78° angle from the seat		O
DKA2540	82° angle from the seat		O
DKA2560	86° angle from the seat		O
DKA2580	90° angle from the seat		•
DKA2590	94° angle from the seat		0
SEAT UPHOLS	TERY*		
DKA2600	Standard Nylon Seat Upholstery with one 5cm central band (Loop)		•
DKA2610	Tension Adjustable Seat Upholstery with one 5cm central band (Loop)		£37
DKA2630	Light Seat Upholstery, mesh material to reduce weight		£57
SEAT CUSHION			
DKA2700	Standard Seat Cushion, 60 mm Double Foam with Visco Top Layer		£55
DKA2720	Matrx Libra Cushion, Black Smooth Cover		£323
DKA2730	Matrx Vi Cushion, Black Smooth Cover - High Risk		£299
SEAT BAG			
DKA2800	Seat Bag, Click in and removable		£61
FOOTREST*			
DKA3200	Footrest Tube, Rigid - Aluminium Hydroformed Profile		*
DKA3210	Footrest Tube, Rigid - Titanium (Rounded tube)		£153
DKA3220	Carbon Footplate, Angle Adjustable	(NC 3290)	£187
DKA3230	Footplate Flip-to Side, Angle Adjustable	(NC 3290)	£119
FOOTREST CO	VER		
DKA3290	Footrest Cover		£31
CLOTHESGUA	RDS/MUDGUARDS*		
DKA3420	Carbon Clothesguards, Flip up, Fold in		*
DKA3440	Carbon Mudguards, Removable	(NC 3490, 3495)	£255
ARMREST			
DKA3480	Tubular Armrest. Height Adjustable, Removable, Non- locking, Swing-away and Comfort Padded, no Clothesguard included		£199
DKA3490	Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included		£136
DKA3495	Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, no Clothesguard included		£136

FRAME COLO	JR*		
DKA4000	Space Black (Standard)		O
DKA4010	Anthracite Grey (Standard)		C
DKA4020	Pure White (Standard)	(NC 4220)	C
DKA4030	Carmine Red (Standard)		O
DKA4040	Azure Blue (Standard)		C
DKA4100	Black Glossy (Special Additional 10 day lead time)	(NC 4190)	£25
DKA4110	Ferrite (Special Additional 10 day lead time)	(NC 4190)	£25
DKA4120	Slate Metallic (Special Additional 10 day lead time)	(NC 4190)	£25
DKA4130	Silver Grey (Special Additional 10 day lead time)	(NC 4190)	£25
DKA4140	Royal Blue (Special Additional 10 day lead time)	(NC 4190)	£25
DKA4150	Ruby Red (Special Additional 10 day lead time)	(NC 4190)	£25
DKA4160	Smokey Grey (Special Additional 10 day lead time)	(NC 4190)	£25
DKA4170	Dynamite Orange (Special Additional 10 day lead time)	(NC 4190)	£25
DKA4180	Forest Green (Special Additional 10 day lead time)	(NC 4190) (NC 4190)	£25
DKA4050	Personal frame colour. Provide RAL/Colour of choice. Price on application	(140 + 150)	POA
FRAME DECAL	_ TYPE*		
DKA4190	Bold Decal (Full coloured transfer on frame)	(NC 4300, 4310)	£25
DKA4195	Tag (Small Küschall brand transfer on front of frame)	(NC 4200, 4210, 4220, 4230, 4240, 4250)	•
BOLD DECAL			
DKA4200	Black		0
DKA4210	Dark Grey		0
DKA4220	Light Grey		0
DKA4230	Red		O
DKA4240	Blue		О
DKA4250	Gold		O
DKA1111	No bold colour, select only with Option DKA4195 Tag	(NC 4190)	O
TAG COLOUR	•		
DKA4300	Black		0
DKA4310	Red		Ο
DKA2222	No tag colour, select only with Option DKA4190 Bold Decal	(NC 4195)	0
FRONT CASTO)RS*		
DKA5000	75mm x 25mm (3" x 1") Supersport		Ο
DKA5050	100mm x 30mm (4" x 1.2") Everyday Black hub and tyre		O
DKA5100	100mm x 30mm (4" x 1.2") Comfort Black hub and tyre		Ο
DKA5150	100mm x 25mm (4" x 1") Starec Black hub and tyre		£69
DKA5200	100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre		£77
DKA5250	100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre		£77
DKA5300	125mm x 30mm (5" x 1.2") Everyday Black hub and tyre		0
DKA5350	125mm x 30mm (5" x 1.2") Comfort Black hub and tyre		O
DKA5400	125mm x 25mm (5" x 1") Starec Black hub and tyre		£69
DKA5450	125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre		£77
DKA5500	125mm x 36mm (5" x 1.4") Aluminium Softroll Silver hub, Black tyre		£77

	POSITION*		
DKA6020	Standard, Fixed Rear Wheel Position	(NC 6100, 6110, 6120, 6130, 6140, 6230, 6240)	•
DKA6030	Light Rear Wheel Extension (3 Positions)	(NC 6070, 6130, 6140, 6230, 6240)	£93
DKA6040	Rear Wheel Extension (5 Positions)	(NC 6070)	£108
CAMBER AXLE	TYPE*		
DKA6050	Vario-Ax, Adjustable Camber can achieve 1°, 3°, 7° and 10°. Only active/scissor brake will lock 10° camber position	(NC 6200, 6210, 6220, 6270, 7230)	£237
DKA6060	Standard Aluminium Axle		•
DKA6070	Carbon Axle		£175
POSITION OF A	AXLE		
DKA6100	1st Position from the front. If 2nd axle is selected, 1st axle must be placed in 1st position and 2nd in the 5th position		0
DKA6110	2nd Position from the front	(NC 6230, 6240)	0
DKA6120	3rd Position from the front	(NC 6230, 6240)	0
DKA6130	4th Position from the front	(NC 6230, 6240)	0
DKA6140	5th Position from the front	(NC 6230, 6240)	0
	LE*		
DKA6200	0° Camber	(NC 6270)	0
DKA6210	3° Camber		•
DKA6220	6° Camber	(NC 6240, 6250)	0
2ND CAMBER	AXLE TYPE. If a 2nd axle is selected, 1st axle must be place	ed in 1st position and the 2nd axle in the 5th position	
DKA6230	Standard Aluminium Axle		£191
DKA6240	Vario-Ax, Adjustable Camber can achieve 1°, 3°, 7° and 10°	(NC 6250, 6260, 6270, 7230)	£237
	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle.		
CAMBER ANG	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle.		
CAMBER ANGI 0° or 3° camber	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle.		oossible in
CAMBER ANGI 0° or 3° camber DKA6250	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle.		oossible in O
CAMBER ANGI 0° or 3° camber DKA6250 DKA6260 DKA6270	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle. 0° Camber 3° Camber		oossible in O
CAMBER ANGI 0° or 3° camber DKA6250 DKA6260 DKA6270	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle. 0° Camber 3° Camber 6° Camber		oossible in O
CAMBER ANGI 0° or 3° camber DKA6250 DKA6260 DKA6270 QUAD QUICK F	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle. 0° Camber 3° Camber 6° Camber RELEASE AXLE Quad Quick Release/Tetra Clip	If only 1 of the axles is a Vario-Axle, 2nd axle is only p	oossible in O O
CAMBER ANGI 0° or 3° camber DKA6250 DKA6260 DKA6270 QUAD QUICK F DKA6300	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle. 0° Camber 3° Camber 6° Camber RELEASE AXLE Quad Quick Release/Tetra Clip	If only 1 of the axles is a Vario-Axle, 2nd axle is only p	oossible in O O
CAMBER ANGI 0° or 3° camber DKA6250 DKA6260 DKA6270 QUAD QUICK F DKA6300 REAR WHEEL	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle. 0° Camber 3° Camber 6° Camber RELEASE AXLE Quad Quick Release/Tetra Clip SIZE*	If only 1 of the axles is a Vario-Axle, 2nd axle is only p	oossible in O O
CAMBER ANGI 0° or 3° camber DKA6250 DKA6260 DKA6270 QUAD QUICK F DKA6300 REAR WHEEL S DKA6410	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle. 0° Camber 3° Camber 6° Camber RELEASE AXLE Quad Quick Release/Tetra Clip SIZE* 24"	If only 1 of the axles is a Vario-Axle, 2nd axle is only p (NC 6550) (NC 6510, 6610, 6640, 6650, 6670, 6680, 6700, 6710, 6720, 6730,	oossible in ○ ○ £71
CAMBER ANGI 0° or 3° camber DKA6250 DKA6260 DKA6270 QUAD QUICK F DKA6300 REAR WHEEL DKA6410 DKA6420	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle. 0° Camber 3° Camber 6° Camber RELEASE AXLE Quad Quick Release/Tetra Clip SIZE* 24" 25" 26"	If only 1 of the axles is a Vario-Axle, 2nd axle is only p (NC 6550) (NC 6510, 6610, 6640, 6650, 6670, 6680, 6700, 6710, 6720, 6730, 6740, 6750, 6760, 6770, 6900, 6910) (NC 0907, 0909, 6520, 6550, 6610, 6630, 6640, 6650, 6670, 6720,	eossible in
CAMBER ANGI 0° or 3° camber DKA6250 DKA6260 DKA6270 QUAD QUICK F DKA6300 REAR WHEEL DKA6410 DKA6420 DKA6430	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle. 0° Camber 3° Camber 6° Camber RELEASE AXLE Quad Quick Release/Tetra Clip SIZE* 24" 25" 26"	If only 1 of the axles is a Vario-Axle, 2nd axle is only p (NC 6550) (NC 6510, 6610, 6640, 6650, 6670, 6680, 6700, 6710, 6720, 6730, 6740, 6750, 6760, 6770, 6900, 6910) (NC 0907, 0909, 6520, 6550, 6610, 6630, 6640, 6650, 6670, 6720,	eossible in
CAMBER ANGI 0° or 3° camber DKA6250 DKA6260 DKA6270 QUAD QUICK F DKA6300 REAR WHEEL DKA6410 DKA6420 DKA6430 REAR WHEEL	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle. 0° Camber 3° Camber 6° Camber RELEASE AXLE Quad Quick Release/Tetra Clip SIZE* 24" 25" 26" TYPE* 24", 26" Küschall Starec Wheel, Silver Cross-Spokes,	If only 1 of the axles is a Vario-Axle, 2nd axle is only p (NC 6550) (NC 6510, 6610, 6640, 6650, 6670, 6680, 6700, 6710, 6720, 6730, 6740, 6750, 6760, 6770, 6900, 6910) (NC 0907, 0909, 6520, 6550, 6610, 6630, 6640, 6650, 6670, 6720,	eossible in
CAMBER ANGI 0° or 3° camber DKA6250 DKA6260 DKA6270 QUAD QUICK F DKA6300 REAR WHEEL DKA6410 DKA6420 DKA6430 REAR WHEEL DKA6510	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle. 0° Camber 3° Camber 6° Camber RELEASE AXLE Quad Quick Release/Tetra Clip SIZE* 24" 25" 26" TYPE* 24", 26" Küschall Starec Wheel, Silver Cross-Spokes, Colourless hub and Black rim	If only 1 of the axles is a Vario-Axle, 2nd axle is only p (NC 6550) (NC 6550) (NC 6510, 6610, 6640, 6650, 6670, 6680, 6700, 6710, 6720, 6730, 6740, 6750, 6760, 6770, 6900, 6910) (NC 0907, 0909, 6520, 6550, 6610, 6630, 6640, 6650, 6670, 6720, 6730, 6740, 6750, 6760, 6760, 6770, 7240)	Possible in
CAMBER ANGI 0° or 3° camber DKA6250 DKA6260 DKA6270 QUAD QUICK F DKA6300 REAR WHEEL DKA6410 DKA6420 DKA6420 DKA6430 REAR WHEEL DKA6510 DKA6520	10° LE 2ND AXLE. Max of 3° difference in camber from 1st axle. 0° Camber 3° Camber 6° Camber RELEASE AXLE Quad Quick Release/Tetra Clip SIZE* 24" 25" 26" TYPE* 24", 26" Küschall Starec Wheel, Silver Cross-Spokes, Colourless hub and Black rim 24", 25" High Performance Wheel Black hub and rim 24", 25" Spinergy LXK, Black hub, rim and spokes in	If only 1 of the axles is a Vario-Axle, 2nd axle is only p (NC 6550) (NC 6510, 6610, 6640, 6650, 6670, 6680, 6700, 6710, 6720, 6730, 6740, 6750, 6760, 6770, 6900, 6910) (NC 0907, 0909, 6520, 6550, 6610, 6630, 6640, 6650, 6670, 6720, 6730, 6740, 6750, 6760, 6760, 6770, 7240) (NC 6650, 6670)	Possible in ○ ○ €71 ◆ ○ ○ €71 ◆ ○ ○ €71

REAR WHEEL 1	YRE*		
DKA6600	24", 25", 26" Schwalbe Rightrun Black/Grey Pneumatic tyre		•
DKA6610	24", 25", 26" Schwalbe Rightrun Blue with K Guard Pneumatic tyre		£40
DKA6620	24", 25", 26" Schwalbe Marathon Plus Smart Guard (Puncture Protected tyre with 5mm layer of flexible rubber)	(NC 6680)	£82
DKA6630	24", 25" Schwalbe One V-Guard Puncture Resistant tyre (Cut-Resistant High-Tech Fibre)		£80
DKA6640	24" Profile Grey Pneumatic tyre 7 Bar		0
DKA6650	24" Solid Profile Grey Tyre	(NC 6680)	£55
DKA6670	24" Solid Light, Grey Tyre	(NC 6680)	£55
ANTI FLAT TYR			
DKA6680	Anti-Flat Tyre Tube Protection		£18
HANDRIM*			
DKA6700	24", 26" Supergripp, mounted at 3 cm distance		£185
DKA6710	24", 26" Steel, mounted at 3 cm distance		£52
DKA6720	24" Titanium, mounted at 3 cm distance		£159
DKA6730	24" Titanium, mounted at 4 cm distance		£159
DKA6740	24" Carbolife Gekko mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip		£265
DKA6750	24" Carbolife Gekko S, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands		£265
DKA6760	24" Carbolife Curve L, mounted at 3 cm distance. Unique shape for low-effort propulsion		£235
DKA6770	24" Carbolife Quadro, mounted at 3 cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim		£295
DKA6780	24", 25", 26" Light Aluminium Grey, mounted at 3 cm distance		•
DKA6790	24", 25", 26" Light Aluminium Black, mounted at 3 cm distance		£34
DKA6800	24", 25", 26" Light Aluminium Grey, mounted at 4 cm distance		0
SPOKEGUARD			
DKA6900	24", 26" Transparent Spokeguard		£104
DKA6910	24", 26" Black Spokeguard		£104
BRAKES*			
DKA7210	Standard Push to Lock Brakes, Lever Bent		•
DKA7230	Performance Brakes, Push to Lock. Brake barrel flips away from tyre	(NC 7300)	£64
DKA7240	Active/Scissor Brakes	(NC 7300)	£64
DKA7250	Light/Compact Active, Aluminium Brakes	(NC 7300)	£172
BRAKE LEVER			
DKA7300	Brake Lever Extension		£41
DKA7500	Anti-Tipper Right		£101
DKA7510	Anti-Tipper Left		£101
DKA7520	Anti-Tipper Both sides		£202
DKA7600	Transportation Kit Includes Posture Belt and Labels	(NC 7800)	£61
DKA7800	Standard Pelvic belt		£61

OTHER OPTIO	NS	
DKA8010	Front Frame Protectors; 2 Piece. Removable with hook and loop	£31
DKA8200	Cane Holder, Mounted Right Side	£62
DKA8210	Cane Holder, Mounted Left Side	£62
DKA8500	Air Pump	£11
DKA8550	Tool Kit	•
DKA8680	24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included	£265

TOTAL PRICE

DKA1310	320 mm
DKA1350	480 mm
DKA1355	500 mm
DKA1410	350 mm
DKA1420	375 mm
DKA1435	400 mm
DKA1475	475 mm
DKA1490	500 mm
DKA1495	525 mm
DKA1500	Position 1 = 140 mm MOST ACTIVE POSITION
DKA1510	Position 2 = 125 mm
DKA1520	Position 3 = 110 mm
DKA1740	450 mm
DKA1745	460 mm
DKA1750	470 mm
DKA1755	480 mm
DKA1760	490 mm
DKA1815	390 mm
DKA1820	400 mm
DKA1825	410 mm
DKA1830	420 mm
DKA1835	430 mm
DKA2050	Titanium
DKA2100	Standard Backrest with Backbrace
DKA2150	Backrest With Lower Backbrace. Backrest Tube Height is 285 mm
DKA2220	Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener
DKA2230	Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener
DKA2240	Without Backrest Upholstery Select if fitting a Matrx/Fixed backrest, no upholstery required
DKA2430	Push Handles Height Adjustable, Rear Set
DKA2440	Push Handles Height Adjustable, Integrated in backrest
DKA2610	Tension Adjustable Seat Upholstery with one 5cm central band (Loop)
DKA2630	Light Seat Upholstery, mesh material to reduce weight
DKA2720	Matrx Libra Cushion, Black Smooth Cover
DKA2730	Matrx Vi Cushion, Black Smooth Cover - High Risk
DKA3210	Footrest Tube, Rigid - Titanium (Rounded tube)
DKA3230	Footplate Flip-to Side, Angle Adjustable
DKA3290	Footrest Cover
DKA3490	Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included
DKA3495	Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, no Clothesguard included
DKA4190	Bold Decal (Full coloured transfer on frame)
DKA4195	Tag (Small Küschall brand transfer on front of frame)
DKA4200	Black
DKA4210	Dark Grey

DKA4220	Light Grov
-	Light Grey
DKA4230	Red
DKA4240	Blue
DKA4250	Gold
DKA4300	Black
DKA4310	Red
DKA5300	125mm x 30mm (5" x 1.2") Everyday Black hub and tyre
DKA5350	125mm x 30mm (5" x 1.2") Comfort Black hub and tyre
DKA5400	125mm x 25mm (5" x 1") Starec Black hub and tyre
DKA5450	125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre
DKA5500	125mm x 36mm (5" x 1.4") Aluminium Softroll Silver hub, Black tyre
DKA6030	Light Rear Wheel Extension (3 Positions)
DKA6040	Rear Wheel Extension (5 Positions)
DKA6070	Carbon Axle
DKA6100	1st Position from the front. If 2nd axle is selected, 1st axle must be placed in 1st position and 2nd in the 5th position
DKA6110	2nd Position from the front
DKA6120	3rd Position from the front
DKA6130	4th Position from the front
DKA6140	5th Position from the front
DKA6200	0° Camber
DKA6210	3° Camber
DKA6220	6° Camber
DKA6230	Standard Aluminium Axle
DKA6240	Vario-Ax, Adjustable Camber can achieve 1°, 3°, 7° and 10°
DKA6250	0° Camber
DKA6260	3° Camber
DKA6270	6° Camber
DKA6420	25"
DKA6430	26"
DKA6510	24", 26" Küschall Starec Wheel, Silver Cross-Spokes, Colourless hub and Black rim
DKA6520	24", 25" High Performance Wheel Black hub and rim
DKA6550	24", 25" Spinergy LXK, Black hub, rim and spokes in Black
SWO0907	24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKüschall
SWO0909	24", 25" Spinergy LX, Chrome hub and 12 spokes in Black. Coloured spokes are available via MyKüschall
DKA6610	24", 25", 26" Schwalbe Rightrun Blue with K Guard Pneumatic tyre
DKA6630	24", 25" Schwalbe One V-Guard Puncture Resistant tyre (Cut-Resistant High-Tech Fibre)
DKA6640	24" Profile Grey Pneumatic tyre 7 Bar
DKA6650	24" Solid Profile Grey Tyre
DKA6670	24" Solid Light, Grey Tyre
DKA6680	Anti-Flat Tyre Tube Protection
DKA6700	24", 26" Supergripp, mounted at 3 cm distance
DKA6710	24", 26" Steel, mounted at 3 cm distance
DKA6720	24" Titanium, mounted at 3 cm distance
DKA6730	24" Titanium, mounted at 4 cm distance

DKA6740	24" Carbolife Gekko mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip
DKA6750	24" Carbolife Gekko S, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands
DKA6760	24" Carbolife Curve L, mounted at 3 cm distance. Unique shape for low-effort propulsion
DKA6770	24" Carbolife Quadro, mounted at 3 cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim
DKA6900	24", 26" Transparent Spokeguard
DKA6910	24", 26" Black Spokeguard
DKA7230	Performance Brakes, Push to Lock. Brake barrel flips away from tyre
DKA7240	Active/Scissor Brakes
DKA7300	Brake Lever Extension
DKA7600	Transportation Kit Includes Posture Belt and Labels
DKA7800	Standard Pelvic belt



Special Instructions Please supply the following information:

Name:

Email:

For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Phone:

Fax:

Description of special request (please include relavent option code or part number):

Customer Ref: Order Number: