




Atom Propel

The Unrivalled Ultra Lightweight



CRASH TEST APPROVED
STANDARD ISO 7176-18

The lightest crash tested self-propel wheelchair available today.



Atom Propel Technical Specifications



Wheelchair Weight

9.1 kg / 20lbs*
Excluding cushions,
anti-tipwheels & footrests



Armrest Height

220mm / 8.6 inches



Overall Size**

Height: 920mm / 36 inches
Length: 1045mm / 42 inches
Width: 635mm / 25 inches



Backrest Height

406mm / 16 inches



Brakes

Attendant – Pull/Push & Lock
User – Push Lock



Folded Flat Size**

Height: 730mm / 29 inches
Length: 1045mm / 42 inches
Width: 270mm / 10 inches



Seat-to-Floor Height

505mm / 20 inches



Weight Capacity

135kg / 21 stone



Tyres - Solid

Front - 152mm / 6 inches
Rear - 559mm / 22 inches



Frame Material

High Strength Aluminium



Seat Width

406mm / 16 inches
457mm / 18 inches
508mm / 20 inches



Footrest

Swing-away
Quick Release



**Crash Test Approved
Standard**
ISO 7176-19:2022

** Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance (+/- 3%). Overall & folded measurements are inclusive of footrests.