



Atom Transit

The Unrivalled Ultra Lightweight



CRASH TEST APPROVED
STANDARD ISO 7176-18

The lightest crash tested transit wheelchair available today.



Atom Transit Technical Specifications



Wheelchair Weight

7.5kg / 16.5lbs*

* Excluding anti-tip wheels, cushions and footrests.



Weight Capacity

135kg / 21 stone



Overall Size**

Height: 920mm / 36 inches
Length: 1030mm / 40 inches
Width: 620mm / 24 inches



Backrest Height

420mm / 16.5 inches



Brakes

Attendant handles
& push to lock



Folded Flat Size**

Height: 730mm / 28 inches
Length: 1030mm / 40 inches
Width: 320mm / 12.5 inches



Armrest Height

210mm / 8 inches



Seat Width

465mm / 18 inches

Seat Depth

400mm / 16 inches



Wheels

Front - 152mm / 6 inches
Rear - 355mm / 14 inches



Seat-to-Floor Height

475mm / 19 inches



Footrest

Swing-away
Quick Release



Crash Test Approved Standard

ISO 7176-19:2022



Frame Material

High Strength Aluminium

** Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance (+/- 3%). Overall & folded measurements are inclusive of footrests.