



Savvy 4 Technical Specifications



Class 2 Mobility Scooter
Can be used on pavements.



Batteries
36Ah x 2



Overall Size
Length: 1285mm / 50 inches
Width: 645mm / 25 inches
Height: 1150mm / 45 inches



Speed Mode
4mph
International speeds available
dependant on jurisdiction and
local legislation.



Braking System
Electromagnetic



Tyres - Pneumatic
Front - 305mm / 12"
Rear - 305mm / 12"



Range
Up to 21 miles*
Range may vary with battery charge/
condition/option, terrain type, tyre
condition & vehicle loading.



Weight Capacity
160kg / 25 stone



Maximum Slope
15 degrees



Motor
Rear Wheel Drive
Sealed Transaxle
500W
24 Volt DC



Scooter Weight
115kg
(with batteries)



Ground Clearance
100mm / 3.9 inches (unloaded)



Turning Circle
1520mm / 60 inches