



Savvy 8 Technical Specifications



Class 3 Mobility Scooter

Can be used on the roads, as well as pavements.



Batteries

45Ah x 2



Overall Size

Length: 1325mm / 52 inches
Width: 652mm / 25 inches
Height: 1220mm / 48 inches



Speed Mode

8mph
International speeds available dependant on jurisdiction and local legislation.



Braking System

Electromagnetic



Tyres - Pneumatic

Front - 305mm / 12"
Rear - 305mm / 12"



Range

Up to 28 miles*
Range may vary with battery charge/condition/option, terrain type, tyre condition & vehicle loading.



Weight Capacity

180kg / 28.3 stone



Maximum Slope

15 degrees



Motor

Rear Wheel Drive
Sealed Transaxle
600W
24 Volt DC



Scooter Weight

125kg
(with batteries)



Ground Clearance

100mm / 3.9 inches (unloaded)



Turning Circle

1530mm / 60.2 inches