

Cabin Car Mk2



Pictured: Limited Edition Cabin Car Mk³

Class	CLASS 3
Weight Capacity	400lbs / 28 Stone 8lbs
Speed¹	8mph (UK & US)*
Ground Clearance^{1,2}	5.25" / 133mm
Length²	56" / 1422mm
Width²	27.25" / 692mm
Height²	56" / 1435mm
Tyres	13" / 330mm
Range¹	75Ah: 22 miles** 100Ah : 30 miles**
Motor	Rear Wheel Drive Sealed transaxle 24 Volt DC Motor 1400w (4 Pole) motor
Brakes	Regenerative Electromechanical
Battery	Dual 12-volt deep cycle 75Ah & 100Ah Options
Battery Charger	8 amp off-board charger
Max Slope²	12 degrees
Seat Back Height²	28.75" / 730mm
Seat Width²	19" / 482mm

* International Speeds:
Available dependant on
jurisdiction and local legislation

**Range is measured on a
smooth surface traveling at
a constant speed. Real world
mileage may be impacted by
terrain and speed.

¹Varies with user weight, terrain
type, battery charge, battery
condition, tyre condition.

²Due to manufacturing
tolerances and continuous
product improvement this
specification can be subject to
variance of +/-3%.

³Limited Edition Cabin Car
pictured. Specification applies
to all Mk 2 variants unless
specified.